

## Conventional Simulation Skin Markers

Identify and differentiate specific points of reference on treatment sites in conventional simulation and on port films with Beekley Medical's SPOT® lead-free metallic pellets, arrows, and line skin markers.

Convenient and easy to use, Beekley Medical's metallic SPOT skin markers for radiation therapy are ideal for marking 3 point set-ups, isocenters, field borders and match lines on treatment sites.

### X-SPOT®

REF 101 Blue QTY 150/box

- 1.5mm lead-free pellet
- easy Pik-UP tab



### Y-SPOT®

REF 102 Green QTY 132/box

- 2.0mm lead-free pellet
- easy Pik-UP tab



### V-SPOT®

REF 603 Florals® QTY 130/box

- 2.5mm lead-free pellet
- easy Pik-UP tab



### D-SPOT®

REF 104 QTY 99/box

- lead-free arrow with easy Pik-UP tab
- for marking nodes and outer canthus



### Z-SPOT®

REF 114 QTY 40/box

- 4.3mm lead-free pellet with easy Pik-UP tab
- for marking port films



### T-SPOT®

REF 111

QTY 275cm cut to measure roll/box

- .8mm diameter metallic line
- for marking scars, nodes, field borders, and match lines



products on this page shown actual size

Call 1-800-233-5539 • Fax 1-800-735-1234 • Visit [www.beekley.com](http://www.beekley.com) • Email [info@beekley.com](mailto:info@beekley.com)

Customers outside the U.S. – Contact your local distributor for pricing and product availability.  
To locate a distributor call +1-860-583-4700 or email [international@beekley.com](mailto:international@beekley.com)